



The diagram shows a light blue circular disc with a dark purple outer border and a white central hole. A red crosshair is centered on the disc, with small red tick marks at the ends of the lines. The disc is divided into four quadrants by the crosshair. Text is located in the top-left quadrant. Below the disc, two horizontal dimension lines are shown. The top one is labeled 'Measures in millimeters' and '≤ 117'. The bottom one is labeled '120'. The disc is enclosed in a thin black rectangular frame.

All bleed areas will be adjusted by our staff based on media which is selected.

Measures in millimeters

≤ 117

120